



**LUNGENE<sup>®</sup>**  
**RAPID TEST**

**COVID-19 IgG/IgM Rapid Test Cassette (WB/S/P)**

A lateral flow chromatographic immunoassay rapid test for the qualitative detection of antibodies (IgG and IgM) to Novel coronavirus in human Whole Blood/Serum/Plasma.

Contains: 25 Test Cassettes      1 Package Insert  
25 Droppers                              1 Buffer

L27VYYMMXXXX      VYYY-MM

IVD   I   CE   Biohazard   4-30°C   40-85°F   30min   15min   5min

**LUNGENE<sup>®</sup>**  
**RAPID TEST**

For in vitro Diagnostic Use Only      Store Sealed at 4-30°C (40-85°F)

IVD   I   CE   Biohazard   4-30°C   40-85°F   30min   15min   5min


COVID-19  
IgG/IgM

ID: \_\_\_\_\_

C  
IgG

IgM

S



**LINGENE COVID-19 IgG/IgM Rapid Test Cassette (WB/S/P) (Rapid)**  
For professional and in vitro diagnostic use only.

**[INTENDED USE]**  
The COVID-19 IgG/IgM Rapid Test Cassette is a lateral flow chromatographic immunoassay for the qualitative detection of antibodies (IgG and IgM) to SARS-CoV-2 in human whole blood/serum/plasma. It provides an aid in the diagnosis of infection with novel coronavirus.

**[SUMMARY]**  
Early January 2020, a novel coronavirus (SARS-CoV-2, formerly known as 2019-nCoV) was identified as the infectious agent causing an outbreak of viral pneumonia in Wuhan, China, where the first cases had their symptoms onset in December 2019. Coronavirus are enveloped RNA viruses that are distributed broadly among humans, other mammals, and birds and that cause respiratory, enteric, hepatic, and neurological diseases. Coronavirus species are known to cause human diseases. Five species—229E, OC-42, NL63, and HKU1—are pathogenic, all of which cause common cold symptoms in humans. Other species, including NL63, HKU1, and 229E, cause severe acute respiratory syndrome (SARS) and Middle East respiratory syndrome (MERS-CoV) and have been shown to cause disease in humans and have been...

**QC PASSED**

**Alcohol Pad**  
70% Isopropyl Alcohol  
For Disinfection Use

**Alcohol Pad**  
70% Isopropyl Alcohol  
For Disinfection Use

**Alcohol Pad**  
70% Isopropyl Alcohol  
For Disinfection Use

**LINGENE RAPID TEST**  
For In Vitro Diagnostic Use Only • Store Sealed at 30°C/86°F

ISO 13485 CE

**COVID-19 IgG/IgM RAPID TEST CASSETTE (WB/S/P)**  
2020030088 EXP: 2022-02

Wangchen Company Limited Co., Ltd.  
162 Youshan Road, Nuhang Subdistrict, Nuhang District, 311221 Hangzhou, China

Wangchen Company Limited Co., Ltd. (Shanghai Branch)  
100000 Shanghai, China



July 30, 2018

Jesse Xia, Manager  
LSI International  
504 E Diamond Ave, Suite 1  
Gaithersburg, MD 20877 US

Re: CR180448

CLIA Waiver for [REDACTED] Test Easy Cup

June 25, 2018

Received: July 5, 2018

CLIA Effective Date: July 30, 2018

### **Categorization Notification (Waived)**

Regulations codified at 42 CFR 493.15 et. seq., implementing the Clinical Laboratory Improvement Amendments of 1988, require the Secretary to provide for the categorization of specific clinical laboratory test systems by the level of complexity. Based upon these regulations, the following commercially marketed test system or assay for the analyte is categorized below:

#### **Test System/Analyte(s): (SEE ATTACHMENT)**

Waived status is applicable to test systems and their instructions approved by the FDA. We recommend that the test system instructions include a statement that the test system is waived under CLIA. Any modification to the test system including test system instructions or a change in the test system name must be submitted to the FDA for the evaluation of waiver. If you change the test system name or your company's name or if a distributor's name replaces your name, you must request another categorization by sending in the revised labeling along with a letter to FDA referencing the document number above.

This complexity categorization is effective as of the date of this notification and will be reported in FDA's CLIA Database: <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfCLIA/search.cfm>. This categorization information may be provided to the user of the commercially marketed test system or assay as specified for the analyte indicated. FDA reserves the right to re-evaluate and re-categorize this test based upon additional information received.

If you have any questions regarding this complexity categorization, please contact Helen Yousefi at [helen.yousefi@fda.hhs.gov](mailto:helen.yousefi@fda.hhs.gov).

Sincerely yours,

Donald St.Pierre  
Acting Director  
Office of *In Vitro* Diagnostics and  
Radiological Health  
Center for Devices and Radiological Health

Regarding the IgG test, we have counted the positive rate of the 77 patients during the convalescence period.

| COVID-19 IgG           |          | Number of patients during the convalescence period | Total |
|------------------------|----------|--|-------|
| CLINICENE <sup>®</sup> | Positive | 75   | 75    |
|                        | Negative | 2  | 2     |
| Total                  |          | 77   | 77    |

#### Cross-Reactivity and Interference

- Other common causative agents of infectious diseases were evaluated for cross reactivity with the test. Some positive specimens of other common infectious diseases were spiked into the Novel coronavirus positive and negative specimens and tested separately. No cross reactivity was observed with specimens from patients infected with HIV, HAV, HBsAg, HCV, IP, H1N1, CMV, FLUA, FLUB, RSV, MP, CP, HPV, etc.
- Potentially cross-reactive endogenous substances including common serum components, such as lipids, hemoglobin, bilirubin, were spiked at high concentrations into the Novel coronavirus positive and negative specimens and tested, separately. No cross reactivity or interference was observed to the device.

| Analytes   | Conc.    | Specimens |          |
|------------|----------|-----------|----------|
|            |          | Positive  | Negative |
| Albumin    | 20mg/ml  | +         | -        |
| Bilirubin  | 20mg/ml  | +         | -        |
| Hemoglobin | 15mg/ml  | +         | -        |
| Glucose    | 20mg/ml  | +         | -        |
| Lipo. Acid | 200μg/ml | +         | -        |
| Lipids     | 20mg/ml  | +         | -        |


- Some other common biological analytes were spiked into the Novel coronavirus positive and negative specimens and tested separately. No significant interference was observed at the lowest level in the table below.

| Analytes             | Conc. (μmol/l) | Specimens |          |
|----------------------|----------------|-----------|----------|
|                      |                | Positive  | Negative |
| Acetaminophen        | 200            | +         | -        |
| Acetic Acid          | 200            | +         | -        |
| Acetylsalicylic Acid | 200            | +         | -        |
| Benzoyl-glycine      | 100            | +         | -        |
| Caffeine             | 200            | +         | -        |
| EDTA                 | 800            | +         | -        |
| Ethanol              | 1.0%           | +         | -        |
| Gentamic Acid        | 200            | +         | -        |
| β-Hydroxybutyrate    | 20,000         | +         | -        |
| Methanol             | 10.0%          | +         | -        |
| Phenolphthalein      | 200            | +         | -        |
| Phenylpropanolamine  | 200            | +         | -        |
| Salicylic Acid       | 200            | +         | -        |

#### Reproducibility

Reproducibility studies were performed for Novel coronavirus-IgG/CLM Rapid Test at three physician office laboratories (POC). Sixty (60) clinical serum specimens, 20 negative, 20 positive and 20 positive, were used in this study. Each specimen was run in triplicate for three days at each POC. The intra-assay agreements were 100%. The inter-site

agreement was 100 %.



Hangzhou Chengene Stech Co., Ltd  
No. 1 Yehuang Road, Yuhang Sub-district, Yuhang District,  
311121 Hangzhou, China

EC REP Shanghai International Holding Corp.GmbH (Europe)  
Elfeldstrasse 80, D-20857 Hamburg, Germany

#### Index of Symbol

|  |   |  |                                     |
|--|---|--|-------------------------------------|
|  | Do not reuse  |  | For in vitro diagnostic use only    |
|  | Store between 4-30°C                                |  | Consult instructions for use        |
|  | Caution   |  | Lot number                          |
|  | Use by  |  | Contains suffragant, for "sp" tests |
|  | Keep away from sunlight                             |  | Keep dry                            |
|  | Manufacturer  |  | Do not use if package is damaged    |
|  | Authorized representative in the European Community |  |                                     |

Version No.: 1.0  
Effective Date: Mar. 04, 2020



**Shanghai International Holding Corporation GmbH (Europe)**

Eiffestrasse 80, 20537 Hamburg Germany

To whom it may concern

**Ihre Zeichen**  
**YourRef.**

**Unsere Zeichen**  
**Our Ref. 02PBJ267**

**Datum**  
**Date: March 9, 2020**

### Declaration of Notification

The undersigned, JIN, Liang, General Manager of Shanghai International Holding Corporation GmbH (Europe), hereby declares that:


  
China

has signed the EC Declaration of Conformity in agreement with the Annex III of the European Directive 98/79/EC on In Vitro Diagnostic Medical Devices and has submitted the required technical documentation, for the following IVD product (for professional use only):

COVID-19 IgG/IgM Rapid Test

The notification to the German Competent Authorities has been carried out on March 9, 2020 by Shanghai International Holding Corporation GmbH (Europe), the appointed Authorized Representative of Hangzhou Clongene Biotech Co., Ltd.

JIN, Liang  
General Manager  
Shanghai International Holding Corporation GmbH (Europe)  
Authorized Representative

  
**Shanghai International Holding Corporation GmbH (Europe)**  
Eiffestrasse 80  
20537 Hamburg

Tel.:(49) 40 2513175  
Fax.:(49) 40 255726  
Mail:  
shholding@hotmail.com

Amtsgericht Hamburg  
HRB 56 583  
Geschäftsführer:  
Liang Jin

Finanzamt Hamburg  
Steuer-Nr.22/795/00590  
Ust-ID-Nr.DE166892350

Bankverbindungen  
SEB AG Hamburg  
BLZ 200 101 11  
Konto:1343865200

# ***EC Declaration of Conformity***

Name and address  
of the manufacturer: [Redacted] Co., Ltd  
Yuhang  
311121 Hangzhou  
China

We declare under our sole responsibility that

the medical device: **COVID-19 IgG/IgM Rapid Test Cassette  
(WB/S/P)-ICOV4212**

of class: **Other**  
according to article 9 of directive 98/79/EC

meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it. The declaration is valid in connection with the "final inspection report" of the device.

Conformity assessment **Directive 98/79/EC Annex III**  
procedure:

Name and address of  
the Authorised Repre-  
sentative **Shanghai International Holding Corporation GmbH  
(Europe)  
Eiffestrasse 80  
20537 Hamburg  
Germany**

[Redacted] **CE**  
[Redacted]  
[Redacted] Representative

Hangzhou, Mar.04.2020  
Place, date

[Redacted]  
Name and function  
*Z. Zhou*